

# HIPOKRAT 2000 REVISION KNEE REPLACEMENT SYSTEM

5 diameters (10, 12, 14, 16, 18)  
in lengths from 90, 120 and  
150 mm. suitable for tibial  
and femoral components  
with diaphyseal slots  
for improved intramedullary  
vascularisation

Offer 9 sizes for  
optimal bone loss coverage

2 thickness option for distal and  
posterior augment provide a  
solution to address bone loss

9 sizes PE-inserts in  
thicknesses from  
10 mm. to 24 mm.

Safety bolt final fixation of the  
stem and tibial metal back

9 sizes metal back  
offers optimal fit

9 sizes and 3 straight (4, 8 and 12 mm.)  
wedge provide an optimal address to bone loss

2, 4, and 6 mm.  
option offset adapters for  
precise component positioning

