




4.5 mm. Düz Rekonstrüksiyon Plak

4.5 mm. Straight Reconstruction Plate



Material: Stainless Steel

Code Kod	Holes Delik	Length Uzunluk
446.0201	2	22 mm
446.0202	3	34 mm
446.0203	4	46 mm
446.0204	5	58 mm
446.0205	6	70 mm
446.0206	7	82 mm
446.0207	8	94 mm
446.0208	9	106 mm
446.0209	10	118 mm
446.0210	11	130 mm
446.0211	12	142 mm
446.0212	13	154 mm
446.0213	14	166 mm
446.0214	15	178 mm
446.0215	16	190 mm
446.0216	17	202 mm
446.0217	18	214 mm

4.5 mm Straight Reconstruction Plate
4.5 mm. Düz Rekonstrüksiyon Plak



Hipokrat Tıbbi Malzemeler İmalat ve Pazarlama A.Ş.

📍 7407/6 Sokak No:10 35060 Pınarbaşı Bornova - İZMİR / TÜRKİYE

☎ + (90) 232 479 11 09

✉ info@hipokrat.com.tr

☎ + (90) 232 479 44 29

🌐 www.hipokrat.com.tr