



Material: Stainless Steel



DCS Plate  
DCS Plak

Code Kod	Holes Delik	Length Uzunluk
497.0101	4	82 mm
497.0102	5	98 mm
497.0103	6	114 mm
497.0104	7	130 mm
497.0105	8	146 mm
497.0106	9	162 mm
497.0107	10	178 mm
497.0108	11	194 mm
497.0109	12	210 mm
497.0110	13	226 mm
497.0111	14	242 mm
497.0112	15	258 mm
497.0113	16	274 mm

Material: Stainless Steel



Lag Screw  
Lag Vidası

Code Kod	Length Uzunluk	Code Kod	Length Uzunluk	Code Kod	Length Uzunluk
496.0121	45 mm	496.0107	80 mm	496.0114	115 mm
496.0101	50 mm	496.0108	85 mm	496.0115	120 mm
496.0102	55 mm	496.0109	90 mm	496.0116	125 mm
496.0103	60 mm	496.0110	95 mm	496.0117	130 mm
496.0104	65 mm	496.0111	100 mm	496.0118	135 mm
496.0105	70 mm	496.0112	105 mm	496.0119	140 mm
496.0106	75 mm	496.0113	110 mm	496.0120	145 mm

Material: Stainless Steel



Compression Screw  
Kompresyon Vidası

Code  
Kod  
496.0300





Hipokrat Tıbbi Malzemeler İmalat ve Pazarlama A.Ş.

📍 7407/6 Sokak No:10 35060 Pınarbaşı Bornova - İZMİR / TÜRKİYE

☎ + (90) 232 479 11 09

✉ info@hipokrat.com.tr

☎ + (90) 232 479 44 29

🌐 www.hipokrat.com.tr