



**Humerus Epikondil Plak**

**Humeral Epicondylar Plate**

Material: Stainless Steel



Humerus Epicondylar Plate  
Humeral Epikondil Plak

Left Sol	Left Sol	Holes Delik	Length Uzunluk
457.0205	457.0105	5	60 mm
457.0206	457.0106	6	80 mm
457.0201	457.0101	7	100 mm
457.0207	457.0107	8	120 mm
457.0202	457.0102	9	140 mm
457.0208	457.0108	10	160 mm
457.0203	457.0103	11	180 mm
457.0204	457.0104	13	200 mm



Hipokrat Tıbbi Malzemeler İmalat ve Pazarlama A.Ş.

📍 7407/6 Sokak No:10 35060 Pınarbaşı Bornova - İZMİR / TÜRKİYE

☎ + (90) 232 479 11 09

✉ info@hipokrat.com.tr

☎ + (90) 232 479 44 29

🌐 www.hipokrat.com.tr