





Screw to be used	: Fixation with Ø2.0 mm. mini cortical screws
Plate profile	: 1.0 mm. x 5.0 mm. with radius
Plate thickness	: 1.0 mm.
Plate width	: 5.0 mm.
Hole spacing	: 6.0 mm.

Material: Stainless Steel



Code Kod	Holes Delik	Length Uzunluk
506.0101	3	17 mm
506.0102	4	23 mm
506.0103	5	29 mm
506.0104	6	35 mm
506.0105	7	41 mm
506.0106	8	47 mm

2.0 mm. Mini Straight Plate (with radius)
2.0 mm. Mini Düz Plak (radyuslu)



Plate profile	: 1.0 mm. x 5.0 mm. (4 to 6 holes)
Plate profile	: 1.5 mm. x 5.0 mm. (7 and 8 holes)
Plate thickness	: 1.0 - 1.5 mm.
Plate width	: 5.0 mm.
Hole spacing	: 5.0 mm.
Screw to be used	: Fixation with Ø2.0 mm. mini cortical screws

Material: Stainless Steel



Code Kod	Holes Delik	Length Uzunluk
506.0601	4	22 mm
506.0602	5	27 mm
506.0603	6	32 mm
506.0604	7	37 mm
506.0605	8	42 mm

2.0 mm. Mini Compression Plate
2.0 mm. Mini Kompresyon Plak



Screw to be used : Fixation with Ø2.0 mm. mini cortical screws
Length 50 mm. can be cut to the required length with wire cutter.
Plate profile : 1.2 mm. x 5.0 mm.
Plate thickness : 1.2 mm.
Plate width : 5.0 mm.
Hole spacing : 5.0 mm.

Material: Stainless Steel



2.0 mm. Mini T-Plate
2.0 mm. Mini T-Plak

Code Kod	Holes in Shaft Şaft Delik	Holes in Head Baş Delik	Length Uzunluk
506.0401	9	3	50 mm.
506.0402	9	4	50 mm.



This plate can bent three-dimensionally using two plates bending forcep also can be cut to the required length with wire cutter.

Plate profile : 0.9 mm. x 5.0 mm.
Plate thickness : 0.9 mm.
Plate width : 5.0 mm.
Hole spacing : 5.0 mm.
Screw to be used : Fixation with Ø2.0 mm. mini cortical screws

Material: Stainless Steel




2.0 mm. Mini Adaptation Plate
2.0 mm. Mini Adaptasyon Plak

Code Kod	Holes Delik	Length Uzunluk
506.0500	20	100 mm



- Screw to be used : Fixation with Ø2.0 mm. mini cortical screws.
Can be cut to the required length with wire cutter.
- Plate profile : 1.2 mm. x 5.0 mm.
- Plate thickness : 1.2 mm.
- Plate width : 5.0 mm.
- Hole spacing : 5.0 mm.

Material: Stainless Steel



Code Kod	Holes Delik	Length Uzunluk
506.0701	20	100 mm

2.0 mm. Mini Straight Plate
2.0 mm. Mini Düz Plak



- Screw to be used : Fixation with Ø1.5 mm. and 2.0 mm. mini cortical screws
- Plate profile : 0.9 mm. x 11.0 mm.
- Plate thickness : 0.9 mm.
- Overall size : 11.0 mm. x 12.0 mm.
- Hole spacing : 7.0 mm. x 6.0 mm.

Material: Stainless Steel



Code Kod
506.0300

2.0 mm. Mini H Plate
2.0 mm. Mini H Plak



- Screw to be used : Fixation with $\varnothing 2.7$ mm. mini cortical screws
- Plate profile : 1.2 mm. x 16.0 mm. with radius
- Plate thickness : 1.2 mm.
- Plate width : 16.0 mm.
- Hole spacing : 12.5 mm.



Material: Stainless Steel

Code Kod	Holes Delik	Length Uzunluk
506.0900	7	45 mm

2.7 mm. Multi Fragment Plate
2.7 mm. Çoklu Kırık Plak



- Screw to be used : Fixation with $\varnothing 2.0$ mm. mini cortical screws

Material: Stainless Steel



2.0 mm. Mini Condylar Plate
2.0 mm. Mini Kondilar Plak

$\varnothing 2.0$ mm.	
Code Kod	Alignment Hizalama
506.1102	left pin
506.1101	right pin



Screw to be used : Fixation with $\varnothing 2.0$ mm. mini cortical screws
Plate thickness : 1.0 mm.
Plate width : 5.0 mm.
Hole spacing : 6.0 mm.



90° Mini L-Plate
90° Mini L-Plak

Left Sol	Right Sağ
472.0301	472.0302



110° Mini L-Plate
110° Mini L-Plak

Left Sol	Right Sağ
472.0101	472.0102



Screw to be used : Fixation with $\varnothing 2.0$ mm. mini cortical screws
Plate thickness : 1.0 mm.
Plate width : 5.0 mm.
Hole spacing : 6.0 mm.

Material: Stainless Steel



Mini T-Plate
Mini T-Plak

Code Kod
472.0200



Hipokrat Tıbbi Malzemeler İmalat ve Pazarlama A.Ş.

📍 7407/6 Sokak No:10 35060 Pınarbaşı Bornova - İZMİR / TÜRKİYE

☎ + (90) 232 479 11 09

✉ info@hipokrat.com.tr

☎ + (90) 232 479 44 29

🌐 www.hipokrat.com.tr