





Slotted to accep tension.

Screw to be used : Fixation with Ø2.7 mm. mini cortical screws
Plate thickness : 2.0 mm.(2 to 6 holes)
2.5 mm.(7 to 12 holes)
Plate width : 8.0 mm.
Hole spacing : 8.0 mm.

Material: Stainless Steel



2.7 mm. Small Fragment Plate
2.7 mm. Küçük Fragment Plak

Code Kod	Holes Delik	Length Uzunluk
470.0101	2	20 mm
470.0102	4	36 mm
470.0103	5	44 mm
470.0104	6	52 mm
470.0105	7	60 mm
470.0106	8	68 mm
470.0107	9	76 mm
470.0108	10	84 mm
470.0109	12	100 mm



Profile one quarter of a 12.0 mm. diameter tube

Plate thickness : 1.0 mm.
Plate width : 7.0 mm.
Hole spacing : 8.0 mm.
Screw to be used : Fixation with Ø2.7 mm. mini cortical screws

Material: Stainless Steel



2.7 mm. Small 1/4 Tubular Plate
2.7 mm. Küçük 1/4 Tübüler Plak

Code Kod	Holes Delik	Length Uzunluk
506.0801	2	18 mm
506.0802	3	26 mm
506.0803	4	34 mm
506.0804	5	42 mm
506.0805	6	50 mm
506.0806	7	58 mm
506.0807	8	66 mm



Screw to be used : Fixation with Ø2.7 mm. mini cortical screws

Material: Stainless Steel



2.7 mm. Small Condylar Plate
2.7 mm. Küçük Kondilar Plak

Ø2.7 mm.	
Code Kod	Alignment Hizalama
506.1002	left pin
506.1001	right pin



Screw to be used : Fixation with Ø2.7 mm. mini cortical screws
Plate Profil : 1.2 mm. x 7.0 mm.
Plate thickness : 1.2 mm.
Plate width : 7.0 mm.
Hole spacing : 8.0 mm.

Material: Stainless Steel



90° Small L-Plate
90° Küçük L-Plak

Left Sol	Right Sağ
471.0101	471.0102



Screw to be used	: Fixation with Ø2.7 mm. mini cortical screws
Plate Profil	: 1.2 mm. x 7.0 mm.
Plate thickness	: 1.2 mm.
Plate width	: 7.0 mm.
Hole spacing	: 8.0 mm.

Material: Stainless Steel



Small T-Plate
Küçük T-Plak

Code

Kod

473.0100



Screw to be used	: Fixation with Ø2.7 mm. mini cortical screws
Plate Profil	: 1.2 mm. x 7.0 mm.
Plate thickness	: 1.2 mm.
Plate width	: 7.0 mm.
Hole spacing	: 8.0 mm.



110° Small L-Plate
110° Küçük L-Plak

Left

Sol

474.0101

Right

Sağ

474.0102



Hipokrat Tıbbi Malzemeler İmalat ve Pazarlama A.Ş.

📍 7407/6 Sokak No:10 35060 Pınarbaşı Bornova - İZMİR / TÜRKİYE

☎ + (90) 232 479 11 09

✉ info@hipokrat.com.tr

☎ + (90) 232 479 44 29

🌐 www.hipokrat.com.tr