



**HRSS Ters Omuz Protezi**

**HRSS Reverse Shoulder Replacement**



without Hole  
Deliksiz

Proximal Humeral Component  
Humerus Üst Gövde

CoCrMo CoCrMo	Titanium Titanyum	Diameter Çap
310.3883	310.3881	ø38 mm.
310.3884	310.3882	ø42 mm.



with Hole  
Delikli

Proximal Humeral Component  
Humerus Üst Gövde

CoCrMo CoCrMo	Titanium Titanyum	Diameter Çap
310.3183	310.3181	ø38 mm.
310.3184	310.3182	ø42 mm.



Humeral Insert  
Humerus Yatak

UHMW-PE UHMW-PE	UHMW-PE (H.C.Linked) UHMW-PE (Y.Ç.Bağlı)	Size Boy
310.3191	310.3201	ø38mm / STD
310.3192	310.3202	ø38mm / +3
310.3193	310.3203	ø38mm / +6
310.3196	310.3206	ø42mm / STD
310.3197	310.3207	ø42mm / +3
310.3198	310.3208	ø42mm / +6



Spacer  
Yükseltici

Code Kod	Size Boy
310.3185	ø38mm / +9mm
310.3186	ø42mm / +9mm



Material: CoCrMo Alloy

Diameter Length Boy	Çap	ø7	ø9	ø11	ø13
90mm.		310.3021	310.3031	310.3041	310.3051
100mm.		310.3022	310.3032	310.3042	310.3052
110mm.		310.3023	310.3033	310.3043	310.3053
120mm.		310.3024	310.3034	310.3044	310.3054
130mm.		310.3025	310.3035	310.3045	310.3055

Cemented Modular Humeral Stem  
Çimentolu Modüler Humerus Stem

Material: Titanium (Tips+HA) Coated

Diameter Length Boy	Çap	ø7	ø9	ø11	ø13
90mm.		310.3101	310.3111	310.3121	310.3131
100mm.		310.3102	310.3112	310.3122	310.3132
110mm.		310.3103	310.3113	310.3123	310.3133
120mm.		310.3104	310.3114	310.3124	310.3134
130mm.		310.3105	310.3115	310.3125	310.3135

Unemented Modular Humeral Stem  
Çimentosuz Modüler Humerus StemBone Plug  
Tıkaç

Diameter Çap	Code Kod
9 mm.	310.3271
11 mm.	310.3272
13 mm.	310.3273
15 mm.	310.3274



Glenosphere  
Glenoid Küresi

Material: CoCrMo Alloy

Standart Standart	Eccentric Eksantrik	Diameter Çap
310.3211	310.3216	ø38 mm.
310.3212	310.3217	ø42 mm.



Metaglène  
Metal Glenoid Yatağı

Material: Titanium Alloy  
Standart Stem

Code Kod	Size Boy
310.3221	ø26mm (for 38mm)
310.3222	ø28mm (for 42mm)



Metaglène  
Metal Glenoid Yatağı

Material: Titanium Alloy  
Long Stem

Code Kod	Size Boy
310.3224	ø26mm (for 38mm)
310.3225	ø28mm (for 42mm)



ø4.5 mm. Locking Metaglène Screw  
ø4.5 mm. Metal Glenoid Yatak Vidası

Material: Titanium Alloy

Cancellous Profile Spongioz Profil	Cortical Profile Kortikal Profil	Size Boy
310.3231	310.3251	16mm
310.3232	310.3252	18mm
310.3233	310.3253	20mm
310.3234	310.3254	22mm
310.3235	310.3255	24mm
310.3236	310.3256	26mm
310.3237	310.3257	28mm
310.3238	310.3258	30mm
310.3239	310.3259	32mm
310.3240	310.3260	34mm
310.3241	310.3261	36mm
310.3242	310.3262	38mm
310.3243	310.3263	40mm
310.3244	310.3264	42mm



**Hipokrat Tıbbi Malzemeler İmalat ve Pazarlama A.Ş.**

📍 7407/6 Sokak No:10 35060 Pınarbaşı Bornova - İZMİR / TÜRKİYE

☎ + (90) 232 479 11 09

✉ info@hipokrat.com.tr

☎ + (90) 232 479 44 29

🌐 www.hipokrat.com.tr