



**Kalkar Destekli Truva Primer Modüler Femur Stem  
Calcar Supported Truva Primary Modular Femoral Stem**



Spout Body  
Kalkar Destekli Gövde

Material : Titanium Alloy  
Cone : 12/14

Sand Blasted Kumlamalı	TIPS TIPS	TIPS+HA TIPS+HA	Size Boy
160.0421	160.0431	160.0441	XSMALL - 47 x 40mm
160.0422	160.0432	160.0442	SMALL - 55 x 45mm
160.0423	160.0433	160.0443	MEDIUM - 60 x 50mm
160.0424	160.0434	160.0444	LARGE - 65 x 55mm

Material : Titanium Alloy



Stem  
Stem

Code Kod	Size Boy
160.0941	ø14 x 115mm
160.0945	ø15 x 115mm
160.0942	ø16 x 115mm
160.0946	ø17 x 115mm
160.0943	ø18 x 115mm
160.0947	ø19 x 115mm
160.0944	ø20 x 115mm
160.0951	ø14 x 150mm
160.0952	ø15 x 150mm
160.0953	ø16 x 150mm
160.0954	ø17 x 150mm
160.0955	ø18 x 150mm
160.0956	ø19 x 150mm
160.0957	ø20 x 150mm
160.0961	ø14 x 190mm
160.0962	ø15 x 190mm
160.0963	ø16 x 190mm
160.0964	ø17 x 190mm
160.0965	ø18 x 190mm
160.0966	ø19 x 190mm
160.0967	ø20 x 190mm

Material : Titanium Alloy



Locking Bolt  
Kilitleme Vidası

Code Kod	Size Boy
160.0451	XSMALL
160.0452	SMALL
160.0453	MEDIUM
160.0454	LARGE





**Hipokrat Tıbbi Malzemeler İmalat ve Pazarlama A.Ş.**

📍 7407/6 Sokak No:10 35060 Pınarbaşı Bornova - İZMİR / TÜRKİYE

☎ + (90) 232 479 11 09

✉ info@hipokrat.com.tr

☎ + (90) 232 479 44 29

🌐 www.hipokrat.com.tr