



Cemented Humeral Component
Çimentolu Humeral Komponent

Material: CoCrMo Alloy

Left Sol	Sağ Right	Diameter Çap
304.1603	304.0103	ø6 x 50 mm.
304.1703	304.0203	ø6 x 70 mm.
304.1803	304.0303	ø6 x 90 mm.
304.1903	304.0403	ø6 x 110 mm.
304.2103	304.0603	ø7 x 50 mm.
304.2203	304.0703	ø7 x 70 mm.
304.2303	304.0803	ø7 x 90 mm.
304.2403	304.0903	ø7 x 110 mm.
304.2603	304.1103	ø8 x 50 mm.
304.2703	304.1203	ø8 x 70 mm.
304.2803	304.1303	ø8 x 90 mm.
304.2903	304.1403	ø8 x 110 mm.



Cemented Ulnar Component
Çimentolu Ulna Komponent

Material: CoCrMo Alloy

Left Sol	Right Sağ
304.1601	304.0101



Boss Core Bar
Orta Mil

Material: Titanium Alloy

Code Kod
303.0104



Plug Screw Axis
Orta Mil Vidası

Material: Titanium Alloy

Code Kod
303.0105



Bush
Bush

Material: UHMWPE

Code Kod
303.0102



Hipokrat Tıbbi Malzemeler İmalat ve Pazarlama A.Ş.

📍 7407/6 Sokak No:10 35060 Pınarbaşı Bornova - İZMİR / TÜRKİYE

☎ + (90) 232 479 11 09

☎ + (90) 232 479 44 29

✉ info@hipokrat.com.tr

🌐 www.hipokrat.com.tr