



Truva Revizyon Modüler Femur Stem
Truva Revision Modular Femoral Stem



Proximal Body
Üst Gövde

Material : Titanium Alloy
Cone : 12/14

STANDART	Sand Blasted Kumlamalı	TIPS TIPS	TIPS+HA TIPS+HA	Size Boy
	160.0401	160.0461	160.0466	50mm
	160.0402	160.0462	160.0467	60mm
	160.0403	160.0463	160.0468	70mm
	160.0404	160.0464	160.0469	80mm

LATERALIZED	Sand Blasted Kumlamalı	TIPS TIPS	TIPS+HA TIPS+HA	Size Boy
	160.0411	160.0471	160.0476	50mm
	160.0412	160.0472	160.0477	60mm
	160.0413	160.0473	160.0478	70mm
	160.0414	160.0474	160.0479	80mm



Stem
Stem

Material : Titanium Alloy

Diameter Çap	190 mm	With Hole	
		225 mm	265 mm
14mm	160.4030	160.4100	160.4700
15mm	160.4040	160.3912	160.3921
16mm	160.4050	160.4200	160.4800
17mm	160.4060	160.3913	160.3922
18mm	160.4070	160.4300	160.4900
19mm	160.4080	160.3914	160.3923
20mm	160.4090	160.4400	160.5000

Material : Titanium Alloy

Locking Bolt
Kilitleme Vidası

Code Kod	Size Boy	Used with Body Length Kullandığı Üst Gövde Boyu
160.5601	M6 x 25mm	50mm
160.5602	M6 x 35mm	60mm
160.5603	M6 x 45mm	70mm
160.5604	M6 x 55mm	80mm

Ø5 mm Fully Threaded
Distal Locking Screw
Ø5 mm Ful Yivli Alt Uç
Kilitleme Vidası

Material : Titanium Alloy

Code Kod	Size Boy
160.5621	Ø5 x 25mm
160.5622	Ø5 x 30mm
160.5623	Ø5 x 35mm
160.5624	Ø5 x 40mm
160.5625	Ø5 x 45mm
160.5626	Ø5 x 50mm
160.5627	Ø5 x 55mm
160.5628	Ø5 x 60mm
160.5629	Ø5 x 65mm
160.5630	Ø5 x 70mm
160.5631	Ø5 x 75mm



Hipokrat Tıbbi Malzemeler İmalat ve Pazarlama A.Ş.

📍 7407/6 Sokak No:10 35060 Pınarbaşı Bornova - İZMİR / TÜRKİYE

☎️ + (90) 232 479 11 09

✉️ info@hipokrat.com.tr

📠 + (90) 232 479 44 29

🌐 www.hipokrat.com.tr