



**HBRS Femur Üst Uç Rezeksiyon Protezi**  
**HBRS Proximal Femur Resection Prosthesis**



Proximal Femur Component  
Femur Üst Gövde

Material: Co. Cr. Mo. Alloy

Code

Kod

648.0101

Material: Titanium Alloy



Puller Bolt  
Çektirme Vidası

Kod

Code

648.0105

Material: Co. Cr. Mo. Alloy



Extension Piece  
Uzatma Borusu

Code Kod	Length Uzunluk
648.0171	40 mm
648.0172	50 mm
648.0173	60 mm
648.0174	70 mm
648.0175	80 mm
648.0176	90 mm
648.0177	100 mm
648.0178	110 mm
648.0179	120 mm

Material: Co. Cr. Mo. Alloy

Cemented Large Conical Stem  
Çimentolu Büyük Konik Stem

Straight Düz	Curved Eğri	Ø
647.1010	647.1001	Ø10 mm.
647.1011	647.1002	Ø11 mm.
647.1012	647.1003	Ø12 mm.
647.1013	647.1004	Ø13 mm.
647.1014	647.1005	Ø14 mm.
647.1015	647.1006	Ø15 mm.
647.1016	647.1007	Ø16 mm.

Material: Titanium Alloy  
TIPS Coated

Straight Düz	Curved Eğri	Ø
647.1030	647.1020	Ø10 mm.
647.1031	647.1021	Ø11 mm.
647.1032	647.1022	Ø12 mm.
647.1033	647.1023	Ø13 mm.
647.1034	647.1024	Ø14 mm.
647.1035	647.1025	Ø15 mm.
647.1036	647.1026	Ø16 mm.

Material: Titanium Alloy  
TIPS+HA CoatedUncemented Hexagonal  
Large Conical Stem  
Çimentosuz Altıgen  
Büyük Konik Stem

Straight Düz	Curved Eğri	Ø
647.1050	647.1040	Ø10 mm.
647.1051	647.1041	Ø11 mm.
647.1052	647.1042	Ø12 mm.
647.1053	647.1043	Ø13 mm.
647.1054	647.1044	Ø14 mm.
647.1055	647.1045	Ø15 mm.
647.1056	647.1046	Ø16 mm.





**Hipokrat Tıbbi Malzemeler İmalat ve Pazarlama A.Ş.**

📍 7407/6 Sokak No:10 35060 Pınarbaşı Bornova - İZMİR / TÜRKİYE

☎ + (90) 232 479 11 09

✉ info@hipokrat.com.tr

☎ + (90) 232 479 44 29

🌐 www.hipokrat.com.tr