



HBRS Ters Omuz Rezeksiyon Protezi
HBRS Reverse Shoulder Resection Replacement



CoCrMo CoCrMo	Titanium Titanyum	Diameter Çap
649.0013	649.0011	ø38 mm.
649.0014	649.0012	ø42 mm.

Proximal Humeral Resection Component
Rezeksiyon Humerus Üst Gövde



Extension Piece ø22
Uzatma Borusu ø22

Material: Co. Cr. Mo. Alloy

Co. Cr. Mo. Co. Cr. Mo.	Titanium Titanyum	Length Uzunluk
649.0071	649.0081	30 mm
649.0072	649.0082	35 mm
649.0073	649.0083	40 mm
649.0074	649.0084	50 mm
649.0075	649.0085	60 mm
649.0076	649.0086	70 mm
649.0077	649.0087	80 mm
649.0078	649.0088	90 mm
649.0079	649.0089	100 mm



Male - Male Adapter
Erkek - Erkek Adaptör

Material: Co. Cr. Mo. Alloy / Titanium Alloy

Co. Cr. Mo. Co. Cr. Mo.	Titanium Titanyum
649.0262	649.0261

**Used for Full Humerus Resection*

**Ful Humerus Rezeksiyon Uygulamalarında*



Humeral Insert
Humerus Yatak

UHMW-PE UHMW-PE	UHMW-PE (H.C.Linked) UHMW-PE (H.C. Linked)	Size Boy
649.0021	649.0031	ø38mm / STD
649.0022	649.0032	ø38mm / +3
649.0023	649.0033	ø38mm / +6
649.0025	649.0035	ø42mm / STD
649.0026	649.0036	ø42mm / +3
649.0027	649.0037	ø42mm / +6



Spacer
Yükseltici

Code Kod	Size Boy
649.0015	ø38mm / +9mm
649.0016	ø42mm / +9mm



Cemented Modular Stem
Çimentolu Modüler Stem

Material: CoCrMo Alloy

Length Boy Çap Diameter	80mm	100mm	120mm
ø7	649.0161	649.0162	649.0163
ø8	649.0165	649.0166	649.0167
ø9	649.0171	649.0172	649.0173
ø10	649.0175	649.0176	649.0177
ø11	649.0181	649.0182	649.0183



Uncemented Modular Stem
Çimentosuz Modüler Stem

Material: Titanium (Tips+HA) Coated

Length Boy Çap Diameter	80mm	100mm	120mm
ø7	649.0185	649.0186	649.0187
ø8	649.0191	649.0192	649.0193
ø9	649.0195	649.0196	649.0197
ø10	649.0201	649.0202	649.0203
ø11	649.0205	649.0206	649.0207
ø12	649.0211	649.0212	649.0213
ø13	649.0215	649.0216	649.0217



Uncemented Splined Modular Stem
Çimentosuz Kanallı Modüler Stem

Material: Titanium Alloy

Length Boy Çap Diameter	80mm	100mm	120mm
ø7	649.0221	649.0222	649.0223
ø8	649.0225	649.0226	649.0227
ø9	649.0231	649.0232	649.0233
ø10	649.0235	649.0236	649.0237
ø11	649.0241	649.0242	649.0243
ø12	649.0245	649.0246	649.0247
ø13	649.0251	649.0252	649.0253



Glenosphere
Glenoid Küresi

Material: CoCrMo Alloy

Standart Standart	Eccentric Eksantrik	Diameter Çap
649.0041	649.0045	ø38 mm.
649.0042	649.0046	ø42 mm.



Metaglene
Metal Glenoid Yatağı

Material: Titanium Alloy
Standart Stem

Code Kod	Size Boy
649.0051	ø26mm (for 38mm)
649.0052	ø28mm (for 42mm)



Metaglene
Metal Glenoid Yatağı

Material: Titanium Alloy
Long Stem

Code Kod	Size Boy
649.0055	ø26mm (for 38mm)
649.0056	ø28mm (for 42mm)



ø4.5 mm. Locking Metaglene Screw
ø4.5 mm. Metal Glenoid Yatak Vidası

Material: Titanium Alloy

Cancellous Profile Spongioz Profil	Cortical Profile Kortikal Profil	Size Boy
310.3231	310.3251	16mm
310.3232	310.3252	18mm
310.3233	310.3253	20mm
310.3234	310.3254	22mm
310.3235	310.3255	24mm
310.3236	310.3256	26mm
310.3237	310.3257	28mm
310.3238	310.3258	30mm
310.3239	310.3259	32mm
310.3240	310.3260	34mm
310.3241	310.3261	36mm
310.3242	310.3262	38mm
310.3243	310.3263	40mm
310.3244	310.3264	42mm



Hipokrat Tıbbi Malzemeler İmalat ve Pazarlama A.Ş.

📍 7407/6 Sokak No:10 35060 Pınarbaşı Bornova - İZMİR / TÜRKİYE

☎️ +(90) 232 479 11 09

✉️ info@hipokrat.com.tr

☎️ +(90) 232 479 44 29

🌐 www.hipokrat.com.tr